

CG1111: Engineering Principles and Practice I

Preparation for Week 7 Studio

Laboratory Oscilloscope Familiarization

-By Dr Henry Tan




Oscilloscope

Oscilloscopes are useful for looking at very fast changes in voltage over time that we could not measure with the multimeter.

■ Channel Buttons* & Signals Colors

- Channel 1 (yellow)
- Channel 2 (blue)
- MATH (pink)

* For toggling the channels/math on & off.




Oscilloscope Probe

- Two Ports

- Positive

- Ground

- ✓ This terminal of both channels are internally connected to ground




- Attenuation Factor

- Ensure 1x is selected for our measurements; **not** 10x, which is meant for signals exceeding the scope's limits, e.g. 50 V

Oscilloscope

■ Test Signal

- For probes compensation and testing
- Approximately 3.2Vpp, 1kHz square wave
- But can also be used to check scope display and settings



Oscilloscope

■ Scale Knobs


○ Vertical scales (volt/div)

✓ 2 mV to 5 V range

✓ Tap the knob to toggle between fine and coarse adjustments

○ Horizontal scales (time/div)


✓ 1 ns to 50 s range



The selected scale settings can be verified from the status bar below the plots.

Oscilloscope

- Position Knobs
 - Vertical position
 - Horizontal position
 - ✓ All three knobs allow quick reset to their respective zero positions by tapping the knob once



Oscilloscope

■ Run Control

○ Auto Set

- ✓ Auto set time/volt scales to best fit of input signal
- ✓ Recommended if signal is unknown


○ Run/Stop

- ✓ Run/stop the waveform

■ Save Display as BMP

○ Print to USB Device

- ✓ USB device has to be connected in advance





Oscilloscope

■ Measure Control



○ Measure button

- ✓ To display the measurements
- ✓ MENU ON/OFF button toggles the menu on & off



○ Entry knob

- ✓ Used together with the softkeys to select a measurement of interest from the submenu
- ✓ Rotate anticlockwise to scroll down the list
- ✓ Tap once to select





Oscilloscope

■ Coupling

○ Channel Coupling

- ✓ **DC coupling** is useful for viewing waveforms that contain AC and DC signals
- ✓ **AC coupling** is useful for viewing waveforms without DC component (i.e. DC offset removed)



- Coupling choices are available in the individual channel menu, invoked by pressing the respective Channel Buttons
- DC coupling is shown with a bar & dots while AC coupling is shown with a tilde sign (~) next to the channel number, located at bottom left corner of the display